

HPLC Analysis of the Semi-Synthetic Penicillin, Amoxicillin

Introduction

Antibiotics are natural substances released by bacteria and fungi which are capable of killing, or inhibiting, competing microbial species. This phenomenon has long been known; but it was not until 1928 that penicillin, the first true antibiotic,

was discovered by Alexander Fleming, Professor of Bacteriology at St. Mary's Hospital in London.¹ However, it 12 years later that pen years later before penicillin was isolated and developed as a medicine by Howard Florey and Ernst Chain.²

Since this breakthrough in therapeutic medicine, derivatives of penicillin were developed with increased efficacy levels. Amoxicillin is one of several such semisynthetic penicillin's. First discovered in the 1960's its patent has now expired. Consequently, amoxicillin and co-amoxiclav is now marketed under many different trade names worldwide. It is commonly used in the treatment of numerous infections including pneumonia, skin infections and Lyme disease.

This application brief illustrates the rapid analysis of amoxicillin, Figure 1, using the Quasar C18 column.

Experimental Conditions

Method Parameters

All LC method parameters are shown in Table 1.

Table 1. LC Method Parameters.

Quasar C18	100 mm	4.6 mm	3 μ m	N9308806
Mobile Phase	A: H ₂ O (+0.1% formic acid) B: ACN (+0.1% formic acid) Gradient: 5-50% B in 10 minutes			
Flow Rate	1 mL/min			
Temp	20 °C			
Wavelength	254 nm			
Injection Vol.	5 μ l			
Analyte	Amoxicillin			

Column Used:

Solvents and Samples

All solvents were HPLC grade and samples were filtered using a 0.45 μ m nylon filter.

Results and Discussion

Amoxicillin, Figure 1, is successfully analysed in just over 5 minutes using the Quasar C18 column, 100 mm in length, Figure 2. Ideally suited to the analysis of small molecules, such as this antibiotic, the Quasar C18 phase provides enough retention via ligand/analyte interactions whilst also maintaining peak shape due to the ultra-high purity silica base and low residual silanol activity.




Figure 1. Chemical Structure of Amoxicillin.




Figure 2. HPLC Analysis of Amoxicillin.

Conclusion

- The Quasar C18 HPLC phase offers high efficiency separation of this semi-synthetic antibiotic.
- The ultra-high purity silica base and low residual silanol activity yields excellent peak shape even for basic analytes, such as amoxicillin.
- Run time could be further reduced by switching to UHPLC and using a Quasar C18 1.7 μ m column.

References

1. Haven KF (1994). *Martell's Handbook of Chemistry and Physics*. 50th Edition. CRC Press, Boca Raton, FL. p. 182. ISBN 978-1-56308-159-0.
2. "Making Penicillin Possible: Norman Heatley Remembers". *ScienceWatch*. Thomson Scientific. 2007. Archived from the original on February 21, 2007. Retrieved 2007-02-13.

Consumables

	Part Number
Nylon filters	02542880

Phase	Length (mm)	I.D. (mm)	µm	Part
Quasar C18	300	3.9	5	N9308800
Quasar C18	250	4.6	5	N9308801
Quasar C18	150	4.6	5	N9308802
Quasar C18	100	4.6	5	N9308803
Quasar C18	50	4.6	5	N9308804
Quasar C18	150	4.6	3	N9308805
Quasar C18	100	4.6	3	N9308806
Quasar C18	50	4.6	3	N9308807
Quasar C18	150	3.0	3	N9308808
Quasar C18	100	3.0	3	N9308809
Quasar C18	50	3.0	3	N9308810
Quasar C18	150	2.1	3	N9308811
Quasar C18	100	2.1	3	N9308812
Quasar C18	50	2.1	3	N9308813
Quasar C18	100	4.6	1.7	N9308814
Quasar C18	50	4.6	1.7	N9308815
Quasar C18	100	3.0	1.7	N9308816
Quasar C18	50	3.0	1.7	N9308817
Quasar C18	100	2.1	1.7	N9308818
Quasar C18	50	2.1	1.7	N9308819

JCANO INGENIERIA S.A. de C.V
 WhatsApp. 55 7129 9832 y Oficina. 55 5444 0793
 Correos: mcانو@jcanoingenieria.com, ventas@jcanoingenieria.com



For a complete listing of our global offices, visit www.perkinelmer.com/ContactUs

Copyright ©2019, PerkinElmer, Inc. All rights reserved. PerkinElmer® is a registered trademark of PerkinElmer, Inc. All other trademarks are the property of their respective owners.